

SNAP TEST PRODUCT INFORMATION



SNAP BLOCK HEATER AND TEST



SNAP TEST FOR ANTIBIOTICS

- 42050 20 TEST SNAP KIT
- 42052 100 TEST SNAP KIT
- 42058 SNAP STARTER KIT
- 42056 SNAP HEATER

The SNAP **Beta-Lactam** test from Idexx Laboratories gets you results in about 10 minutes and tests for Penicillin G, amoxicillin, ampicillin, ceftiofur and cephalixin.

Rely on SNAP tests for accurate, efficient residue testing, all on a single platform.

BETA-LACTAM TEST KIT
Detects penicillin G, amoxicillin, ampicillin, ceftiofur and cephalixin residues in raw, commingled cow milk at or below the U.S. FDA-established tolerance/safe levels. Beta-lactams are the most widely used antibiotics in the treatment of mastitis and other infections in dairy cattle. Some common treatments that contain beta-lactam drugs are ToDAY®, ToMORROW®, Albacillin®, Cefa-Lak®, Quartermaster®, Cefa-Dri®, Hetacin-K® and NAXCEL®.

TETRACYCLINE TEST KIT
Detects tetracycline, chlortetracycline and oxytetracycline residues in raw, commingled cow milk at or below the U.S. FDA-established tolerance/safe levels. Tetracycline antibiotics are used in the treatment of mastitis and other infections in dairy cattle. Some of the most common tetracycline treatments include LA-2000®, Oxytet™ and Maxim-200®.

AFLATOXIN M1 TEST KIT
For the screening detection of aflatoxin M1 residues in raw, commingled cow milk.

GENTAMICIN TEST KIT
For the screening detection of gentamicin residues in raw, commingled cow milk.

SULFAMETHAZINE TEST KIT
For the screening detection of sulfamethazine residues in raw, commingled cow milk.

SNAP® test kits are backed by a dedicated team of technical service professionals available to answer your product, regulatory and industry questions. Call IDEXX Technical Services at 1-800-321-0207.

*See the package insert that accompanies your kit.

IDEXX LABORATORIES
1-800-120-0000

© 2007 IDEXX Laboratories, Inc. All rights reserved. • 7/06-02