

# DELVOTEST PRODUCT INFORMATION



## DELVOTEST P SP/NT SP

Standard diffusion test for the determination of antibiotic residues in milk.

**42010 DELVOTEST P** – 2 ½ hour test for antibiotics. Components of the test kit: 25 ampules containing bacillus stearothermophilus var. calidolactis in solid medium, 25 tablets containing nutrients based on glucose and indicator, dosing syringe for sampling milk, 25 disposable pipettes, tweezers, and complete instructions.

**42043 DELVOTEST SP/NT** – 3 hour test for antibiotics. Tests for the same antibiotics as the Standard “P” Kit plus sulfa compound drugs. Components same as above.

**42020 DELVOTEST STARTER KIT**

**42014** 8 well – 120 volt incubator

**42013** DELVOTEST thermometer

**42011** DELVOTEST SP/NT 100

The Delvo P detects for:

Penicillin(G)  
Cloxacillin  
Ampicillin  
Amoxicillin  
Cephapirin  
Ceftiofur  
Tetracycline  
Oxytetracycline  
Spiramycin  
Gentamicin  
Pirlimycin

And many more...

The SP/NT detects the above, plus the following:

Sulfamethazine  
Sulfadimethoxine  
Sulfathiozole  
Sulfadiazine  
Dapsone



**DELVO P MINI – 25 TESTS**



**SP/NT 100 Kit**